

VITÁLE®

Hydrocolloid Dressings (Odor Eliminating Formula)



Vitále® Hydrocolloid Dressings are now odor-eliminating! Special non-toxic formula encapsulates odor molecules and holds them safely within the dressing. At dressing change time, you throw away the odor molecules too.

Vitále Hydrocolloid (standard thickness) dressings have a special tapered edge that makes the dressing much more comfortable and helps to eliminate roll-up at the edges.

Vitále Thin Hydrocolloid Dressings have the thin profile of other "extra-thin" hydrocolloids, but the superior absorbency of our formulation allows our dressings to be used even for wounds with moderate drainage.

Indications: Hydrocolloids are indicated for partial & full thickness wound (Stages II-IV), with light to moderate exudation. Hydrocolloids are for single use only and can be used for treatment of the following wounds:



Light to Moderate Exuding



Superficial or Deep

- Pressure Ulcers
- Arterial Ulcers
- Post-Operative Wounds
- 1st & 2nd Degree Burns (not appropriate for 3rd Degree)
- Venous Ulcers
- Diabetic Ulcers
- Skin Donor Sites

Features & Benefits:

Odor Eliminating
Absorbent

Holds exudates and their odors deep within the dressing.
Absorbs exudates and holds moisture within the dressing, allowing usage for up to 7 days.

Moist wound environment
Conforms to wound contours

Keeps the wound moist for the ideal healing environment.
Body heat makes the dressing more conformable, and it may be cut to fit.

Sacral Size available

New Sacral Hydrocolloid is convenient for use in sacral wounds (picture to right).

Tapered Edge

The new tapered edge prevents roll-ups at the edges and is much more comfortable.



Product #:	Description:	Qty/box:	HCPC Codes:
21044	Vitále Thin Hydrocolloid Dressing, 4" x 4" (10cm x 10cm)	10	A6234
21222	Vitále Hydrocolloid (standard) Dressing, 2" x 2" (5cm x 5cm)	10	A6234
21244	Vitále Hydrocolloid (standard) Dressing, 4" x 4" (10cm x 10cm)	10	A6234
21290	Vitale Sacral Hydrocolloid Dressing, 6" x 7" (15cm x 18cm)	5	A6235

Information:

Hydrocolloids were first used in products to adhere stoma bags to patients' skin...and clinicians noted that the edges of the stoma under the hydrocolloid looked healthier and showed signs of healing. This led to the development of hydrocolloid dressings for wound care. Properly made hydrocolloid dressings absorb and hold fluids...thereby keeping wounds moist. Vitále Hydrocolloids promote re-hydration and autolytic debridement of wounds, while not disintegrating into the wound bed. Hydrocolloids work best on wounds with light to moderate drainage. Vitále Hydrocolloid Dressings are sterilized by exposure to gamma radiation and must not be re-sterilized.

