

# Glucertain™ Blood Glucose Monitoring System



The **Glucertain™ Blood Glucose Monitoring System** is the state-of-the-art for monitoring insulin therapy for diabetic patients in the Long-Term Care setting.



No Code...you don't have to enter a code into the meter to utilize Glucertain Test Strips. All Glucertain Strips are calibrated to the Glucertain Meters.



GDH Technology...a unique Glucose Dehydrogenase chemistry system eliminates inaccurate readings because of interference from patient oxygen therapy.



5 Second Result...the meter counts down for 5 seconds and displays the reading. Quick and accurate.



Mini Sample...only 0.6 microliters required to obtain a reading. That's a very small droplet of fingertip blood.

Glucertain Blood Glucose Strips, though they incorporate state-of-the-art technology, are available at economical pricing. So you get low price & the latest technology.

**Glucertain Clinical Safety Lancets** are the perfect choice for glucose finger-sticks when safety, comfort & cost are the criteria. Available in 2 sizes: 28 gauge (standard) and 26 gauge for fingertips with thicker skin.

**Safety:** Self-actuated, automatically retractable, use once & dispose safely.

**Comfort:** The small size (28 gauge is standard) means the ultimate in comfort for patients.

**Cost:** Glucertain contracts (on Lancets & Blood Glucose Testing Meters & Strips) are available to extend the lowest cost to facilities.

For more information on Glucertain Blood Glucose Products, call us at 800-892-2391 or email [glucertain@cellerallc.com](mailto:glucertain@cellerallc.com).

## Ordering Data:

### Product #: Description:

30100	Glucertain Blood Glucose Meter
30150	Glucertain Blood Glucose Strips, 25/vial, 2 vials/box
30180	Glucertain Blood Glucose Control Solution, Low
30190	Glucertain Blood Glucose Control Solution, High
30226	Glucertain Safety Lancets, 26 ga, 100/box, 20 boxes/case
30228	Glucertain Safety Lancets, 28 ga, 100/box, 20 boxes/case

