

# VITÁLE®

## Wound Care Accessories



**Vitale® Skin Protectant Wipes** are "No-Sting". They do not contain alcohol and thus do not sting when applied to delicate peri-wound tissue. Skin Protectant Wipes are used to provide a film coating to tender tissue before applying an adhesive dressing. They protect the skin, provide secure adhesion, and make removal easier at dressing change time. They are also often used to add protection to skin subject to abrasion, such as heels, elbows, etc.

### Features & Benefits:

- Formula contains no alcohol . . . . . "No-Sting" wipes prevent patient discomfort when applied.
- Extra-moist wipes . . . . . Apply a thick film barrier for better protection.
- Economical . . . . . Compete favorably with inexpensive alcohol-containing wipes that "sting".

**Vitale® Wound Cleanser/Spray** is a gentle, non-toxic wound cleanser spray that is also an antiseptic. The surfactant used in the Cleanser is a low-sudsing derivative of sugar, and is blended in normal saline, so it is tissue-friendly. Benzalkonium Chloride is an effective anti-microbial, so the spray helps reduce bioburden in the wound. Finally, Vitale Wound Cleanser/Spray contains zinc to help in the healing process. The Spray is designed to leave the wound bed moist when cleansing is complete.



### Features & Benefits:

- Antiseptic formula . . . . . Reduces microbial contamination in the wound bed.
- Low sudsing normal saline formula . . . Gentle action to cleanse the wound thoroughly and safely.
- Contains zinc . . . . . Conditions the wound bed for healing.
- Ergonomic package . . . . . Easy to use & store.



**Vitale® Composite & Cover Dressings** are a primary dressing for minor wounds or a secondary dressing for gels & ointments, calcium alginate dressings, wound fillers, non-adhesive foams, collagen dressings, etc. Cover Dressings are a one-step replacement for gauze & adhesive tape, and provide a soft, conformable surface to protect the wound and tender surrounding tissue. Vitale Composite Dressings add a layer of breathable film on the outside to provide waterproof protection. Vitale Composite & Cover Dressings are sterile and packed in individual peel packages to help prevent cross-contamination.

### Features & Benefits:

- Soft non-woven fabric . . . . . Provide a soft, conformable protective cover.
- Replaces gauze & tape . . . . . Convenient and secure cover for gels and dressings
- Sterile, pouch packaged . . . . . To help prevent contamination of the primary dressing or the wound.

Product #:	Description:	Qty/box:	HCPC Codes:
24004	Vitale Composite Dressings, 4"x4"	30	A6203
24007	Vitale Composite Dressings, 6"x7"	30	A6203
24502	Vitale Cover Dressings, 2"x2"	30	A6219
24504	Vitale Cover Dressings, 4"x4"	30	A6219
24506	Vitale Cover Dressings, 6"x6"	30	A6220
24544	Vitale Cover Dressings, 4"x14"	20	A6221
24950	Vitale Skin Protectant Wipes	100	A5120
25500	Vitale Wound Cleanser Spray, 240 ml	12	A6260

