

VITÁLE®

Calcium Alginate Dressings

Silver Calcium Alginate Dressings



New...improved Vitále® Calcium Alginate Dressings have enhanced gelling properties. Our exclusive new process also provides better fabric uniformity and conformability in both plain and silver.

Calcium Alginate creates an ideal wound healing environment and is indicated for the treatment of wounds with moderate to heavy exudate which may also be infected. Superior absorption extends wear time, and dressing changes are easily completed without disturbing granulating tissue.

Vitále Calcium Alginate Dressings are ideal for draining wounds, as they reduce the risk of maceration to surrounding skin. The flexibility of the dressings & ropes facilitates use in cavity or tunneled wounds.

Vitále Silver Calcium Alginate Dressings and Ropes are ideal for draining wounds that are infected. The flexibility of the dressings & ropes facilitates use in cavity or tunneled wounds.

Indications: Vitále Calcium Alginate Dressings and Ropes with Silver are for single use only and can be used for the treatment of the following wounds:



Moderate to Heavy Exuding



Superficial or Deep

- Pressure Ulcers
- Arterial Ulcers
- Post-Operative Wounds
- 2nd Degree Burns (not indicated for dry wound conditions or for 3rd Degree Burns)
- Venous Ulcers
- Diabetic Ulcers
- Skin Donor Sites

– Ordering Information –

Cat. #:	Description:	Size:	Qty/box:	HCP Code:
20500	Vitále Silver Calcium Alginate Ropes		5	A6199
20522	Vitále Silver Calcium Alginate Dressings	2"x 2" (5cm x 5cm)	10	A6196
20544	Vitále Silver Calcium Alginate Dressings	4"x 4" (10cm x 10cm)	10	A6196
20545	Vitále Silver Calcium Alginate Dressings	4"x 5" (10cm x 12cm)	10	A6197
20548	Vitále Silver Calcium Alginate Dressings	4"x 8" (10cm x 20cm)	5	A6197
20600	Vitále Calcium Alginate Ropes		5	A6199
20622	Vitále Calcium Alginate Dressings	2"x 2" (5cm x 5cm)	10	A6196
20644	Vitále Calcium Alginate Dressings	4"x 4" (10cm x 10cm)	10	A6196
20645	Vitále Calcium Alginate Dressings	4"x 5" (10cm x 12.5cm)	10	A6197

Features & Benefits:

- Anti-Bacterial Silver The Calcium Alginate fibers are impregnated with silver, providing antimicrobial activity to help prevent sepsis.
- Highly absorbent Absorbs many times its own weight, allowing extended wear time depending on the volume of exudate.
- Moist wound environment The action of exudate on alginate fibers produces a gel that keeps the wound moist.
- Facilitates debridement The moist gel helps autolytic debridement.
- Easy to remove Fast and efficient wound clean-up without effect on granulating tissue.
- Conforms to wound contours Soft and conformable, and may be molded or cut to fit wound. As the fibers gel, they make even more intimate contact with the wound bed.

